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Druckversion: 2020-05

Content	Page
1. Safety Instructions	9
1.1 Classification of the Safety Instructions	9
1.2 All Instructions for Your Safety	9
2. Use	10
2.1 Intended Use	10
2.2 Indication	10
2.3 Qualification	11
2.4 Application	11
2.5 Product Range	11
3. Maintenance of the Orthosis	11
3.1 Dirt Removal from the System Joint	11
4. Disposal of the Orthosis	11
5. Signs and Symbols	12
6. CE Conformity	12
7. Legal Information	12

Instructions for Use for Patients

Waterproof Free Moving System Knee Joints

Dear Patient,

You have received an individually produced orthosis with a high quality FIOR & GENTZ system knee joint from your orthotist or a qualified/trained expert.

1. Safety Instructions

1.1 Classification of the Safety Instructions

 DANGER	Important information about a possible dangerous situation which, if not avoided, leads to death or irreversible injuries.
 WARNING	Important information about a possible dangerous situation which, if not avoided, leads to reversible injuries that need medical treatment.
 CAUTION	Important information about a possible dangerous situation which, if not avoided, leads to light injuries that do not need medical treatment.
<i>NOTICE</i>	Important information about a possible situation which, if not avoided, leads to damage of the product.

All serious incidents according to Regulation (EU) 2017/745 which are related to the product have to be reported to the manufacturer and to the competent authority of the Member State in which the orthotist or qualified/trained expert and/or the patient is established.

1.2 All Instructions for Your Safety

DANGER

Potential Traffic Accident Due to Limited Driving Ability

Gather information about all issues concerning safety and security and potential dangers before driving a motor vehicle with orthosis.

WARNING

Risk of Falling Due to Loosened Screws

Do not secure screws for the system joint on your own. Consult your orthotist or a qualified/trained expert about the maintenance intervals to be respected.

WARNING

Risk of Falling Due to Improper Handling

Have your orthotist or a qualified/trained expert inform you about the correct use of the system joint and potential dangers.

WARNING

Risk of Falling Due to Improper Shoe/Wrong Shoe Pitch

Wear a shoe to which your orthosis is adjusted in order to avoid joint dysfunction.

WARNING

Risk of Falling Due to Orthosis Adjustments

If you notice any changes on the orthosis (e.g. loosely attached joint components or changed spring forces), immediately contact your orthotist or a qualified/trained expert.

NOTICE

Limitation of the Joint Function Due to Improper Dirt Removal

Remove dirt from the orthosis and the system joint as described in these instructions for use. Do not grease the system joint on your own. If necessary, consult your orthotist or a qualified/trained expert.

Any modification of the system joint from your side is not permitted.

2. Use

2.1 Intended Use

The FIOR & GENTZ waterproof free moving system knee joints are exclusively for use for orthotic fittings of the lower extremity. The system joint provides motion control and is only allowed to be used for producing a KAFO. Every system joint influences the orthosis' function and thus also the function of the leg. Depending on the type of orthotic fitting and the patient's individual physical abilities, driving a motor vehicle safely should be guaranteed.

2.2 Indication

The indication for the treatment with an orthosis for the lower extremity is a pathological gait. This can be caused, for example, by central, peripheral, spinal or neuromuscular paralyses, structurally conditioned deformities/malfunctions or surgery.

The physical conditions of the patient such as muscle strength or activity level are decisive for the orthotic treatment. A safe handling of the orthosis must be ensured.

2.3 Qualification

The system joint must be handled by an orthotist or qualified/trained expert.

2.4 Application

All FIOR & GENTZ system joints were developed for everyday life activities such as standing and walking. Extreme loads connected to activities like running, climbing and parachuting are excluded. The NEURO CLASSIC H₂O system knee joint is waterproof and, therefore, it is suited for a usage in wet areas.

2.5 Product Range

The following waterproof free moving system knee joints are part of the FIOR & GENTZ product range:

NEURO CLASSIC H₂O

3. Maintenance of the Orthosis

Ask your orthotist or a qualified/trained expert to check the system joint of your orthosis every **6 months**. Bring these instructions for use to each follow-up and let your orthotist or a qualified/trained expert enter the next appointment in the table on the back. Never carry out maintenance work or other adjustments and repairs yourself. In the case of children and people with cognitive impairments, we would like to point out to you as parents or care team that you must regularly check the orthosis and the system joint for signs of wear. If you notice any changes, immediately contact your orthotist or a qualified/trained expert.

3.1 Dirt Removal from the System Joint

Remove dirt from the system joints on a regular basis. Use a dry cloth and clean the system joint only superficially. Remove visible dust and lint from the mechanics by using tweezers. Check the orthosis in straight and flexed position.

In order to optimise the lifespan, we recommend rinsing the orthosis with water, especially after using it in salt water, chlorine water and sand.

4. Disposal of the Orthosis

If you no longer need the orthosis, please return it to your orthotist or a qualified/trained expert. The product must not be disposed of with the residual waste (fig. 1).


fig. 1

5. Signs and Symbols

Symbols on the Packaging

medical device

6. CE Conformity

We declare that our medical devices as well as our accessories for medical devices are in conformity with the requirements of Regulation (EU) 2017/745. Therefore, the FIOR & GENTZ products bear the CE marking.

7. Legal Information

With the purchase of this product, our General Terms and Conditions of Business Transactions, Sales, Delivery and Payment will apply. The warranty expires, for example, if the product is mounted several times. Please note that the product is not supposed to be combined with other components or materials than with those recommended by the FIOR & GENTZ Orthosis Configurator. Should the product be combined with products from other manufacturers, a written consent by the seller is required.

The information in these instructions for use is valid at the date of printing. The contained product information serves as a guideline. Subject to technical modifications.

All rights, particularly the distribution, copy and translation of this manual or any part of it, in paper or as electronic document, must be authorised in writing by FIOR & GENTZ Gesellschaft für Entwicklung und Vertrieb von orthopädietechnischen Systemen mbH. Reprints, copies and any other electronic reproduction, even partial, must be authorised in writing by FIOR & GENTZ Gesellschaft für Entwicklung und Vertrieb von orthopädietechnischen Systemen mbH.

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